

# IEEE NPSS Real Time School

# Time-of-Flight Data Acquisition Lab

Stefan Ritt Aziza Zendour June 2024

SSID:	TOF
PW:	IEEENPSS
URL:	http://pi .local



# 1 Introduction

This lab teaches you how to measure the speed and direction of cosmic muons generated in the upper atmosphere of the earth with the Time-of-flight (ToF) technique. Two scintillator plates are read out by silicon photomultipliers (SiPM) each on two sides. In a first measurement, the time difference between the signals on both plates is measured to determine the speed of the muons, while in a second measurement the two plates are rotated to measure the direction of the muons.



# 2 Equipment

# 2.1 Silicon Photomultiplier (SiPM)

SiPM are photon-counting devices made up of multiple avalanche photodiodes (APD) working in Geiger mode. A photon impinging on the silicon surface creates an electron-hole pair through ionization. These carries are then accelerated by a high voltage and create secondary electron-hole pairs (avalanche effect). This process continues and give a gain of 10<sup>5</sup>-10<sup>6</sup>, so a single photon can generate a measurable electrical signal.



Figure 1: Working principle of an avalanche photodiode (source: Hamamatsu MPPC Data Sheet)

The APD is operated in Geiger mode, so the device goes in breakdown mode where it becomes conductive even if operated in reverse polarity. To prevent damage of the diode, a quench resistor is used to in series (see Figure 2). If the current goes above a certain limit, the voltage drop across the resistor becomes so big that the voltage across the diode drop below the breakdown voltage and the avalanche stops. After a certain recovery time, the diode becomes sensitive again. A single diode goes into breakdown independent from the number of primary photons and is therefore not suited to measure the amount of light and therefor the energy a particle deposits in a scintillator. To overcome this problem, several APDs are combined on a single chip. This arrangement is called Silicon Photomultiplier (SiPM) or Multi-Pixel Photo Detector (MPPC).



Figure 2: SiPM with many APDs each having its own quench resistor. The right picture shows a photograph of a real SiPM.

# 2.2 Scintillator



#### Figure 3: Scintillator coupled to two SiPMs

## 2.3 WaveDREAM Board

The WaveDREAM board is an integrated data acquisition board developed at the Paul Scherrer Institute, Switzerland. It contains 2 x DRS4 chips which sample the input signal with up to 5 GS/s. It has amplifiers with variable gain for each channel to digitize directly signals from silicon photomultipliers (SiPM). An integrated high voltage board devlivers the ~50 V bias voltage for the SiPMs. The board is read out via Gigabit Ethernet to a PC or a RaspberryPi computer.



Figure 4: Schematics of the WaveDREAM Board with 16 input channels, ADC, DRS4 and FPGA



Figure 5: [Left side]: WaveDREAM board (bottom) with 240V high voltage Cockroft-Walton Generator piggy-back (top). [Right side]: Boxed WaveDREAM boards

## 2.4 RaspberryPi



Figure 6: PaspberryPi Computer

The RaspberryPi Computer runs the WaveDREAM Server Program (wds) to receive the individual data packets from the WaveDREAM board and make the waveforms available to a web browser based oscilloscope application.

## 2.5 WDS Oscilloscope



The oscilloscope application has built-in measurement functions with which one can analyze the signals coming from particle detectors.

## 2.6 I-V Curve

The current flowing through a SiPM depends on the voltage. The higher the voltage, the more ionization inside the silicon and the higher the current. At a certain voltage, the Geiger-Mueller

region starts, where we have full ionization of a SiPM pixel. This is called the "breakdown-voltage"  $V_b$ . One can measure this breakdown voltage by recording a curve of current vs. voltage (see Figure 7). If time is left, the I-V curve can be recorded by changing the high voltage and reading the current for each channel as described later. The typical operation voltage is the  $V_{op} = V_b + 3V$ . Please make sure not to exceed a current of 1.5 uA!



Figure 7: Typical I-V curve of a SiPM

# 3 Measurement

## 3.1 Preparation

Please make sure all cables and devices are connected as shown in the following picture. The lower scintillator plate should be connected to channel 0 and 1 (CH00, CH01) to the silver data acquisition board, and the upper plate must be connected to CH02 and CH03.

**Note:** Channels 0 and 1 (CH00, CH01) can sometimes be noisy. To reduce noise, you may switch the connections as follows: CH00 -> CH02, CH01 -> CH03, CH02 -> CH04, and CH03 -> CH05.



Now connect your browser to the RaspberryPi computer of your setup using the URL printed on the cover page of this manual. Make sure not to connect to the Pi of your neighbour! Then, configure the DAQ board according to following scheme (note the blue labels on the pictures below):

- 1. Turn channel 0, 1, 2, 3 on, all other channels off (must be grey)
- 2. Click on "Channels" to reveal the channel configuration panel. Set the HV of the first four channels to the value printed on the detectors. Usually either 54 V or 56 V. Make sure to have some current flow of 0.2-0.7  $\mu$ A. If there is no current, report to your instructor.
- 3. Click on "Config" to reveal the configuration panel
- 4. Set the Trigger Level to -45 mV and the Delay to 94 ns
- 5. Set the Trigger Type to "normal"
- 6. Set the Gain to 50



•		•			Cha	nnel Co	nfiguration				
		Chn	Gain	PZC	Trigger	Level	HV		Curren	t	
		0	50 ~		-32	mV	54	V	0.506	uA	
ľ		1	50 ~			nV	54	V	0.314	uA	
		2	50 ~			nV	54	V	0.628	uA	
		3	50 ~		-32	mV	54	V	0.419	uA	
		4	50 ~		-32	mV	0	V	0	uA	
		5	50 ~		-32	mV	0	V	0	uA	
		6	50 ~		-32	mV	0	V	0	uA	
		7	50 ~		-32	mV	0	v	0	uA	
1		8	50 ~		-32	mV	0	V	0	uA	~~~
	-	9	50 V		-32	mV	0	V	0	υA	~~~~

Next configure the trigger. Click on "Trigger Pattern" to reveal the trigger configuration page blow.

~~	~~~	~~~	~~~~			~	+	~		1	-	~		~~							
	Trigger Patterns																				
<u> </u>			_	~					lear		Or	All		Help	•						
~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	Chn	Pol	P00	Р	1 P	02	P03	P04	P05	P06	P07	P08	P09	P10	P11	P12	P13	P14	P15	P16	P17
	СНО	-	٠																		
	CH1	-	٠	Ц	4																
	CH2	-	٠	Ц	_	_															
τ	СНЗ	-	•	J	+																
	CH4	+	_	γ	_	_															
	CH5	+																			
	CH6	+			4	_		_								_					
	CH7	+																			
	CH8	+				_															
	CH9	+																			
	CH10	+																			
	CH11	+																			
	CH12	+			_	_															
	CH13	+																			
	CH14	+		L	+	_		_								_					
	CH15	+			4																
	EXT	+																			

Here, switch the polarity of CH0-CH3 to negative ("-") and turn the coincidence on (solid dot under the "P00" header). This requires all four channels to carry a signal to trigger the board. Make sure all other Pxx columns are grey (inactive).

Now define the time measurement. Click on "Measure", then click on "+" to add a new measurement. Select "Chn delay", then "CHO", Thr1: -20, "CH2", Thr2: -20. This defines a channel delay measurement between channels 0 and 2 (not 1!) which are the lower and upper scintillator plates. The threshold for the time measurement is set to -20 mV. Turn on "Display Histograms". If everything works correctly, you should see a peak around "0" in the bottom histogram after a few minutes. Now use the "handles" left and right of the histogram to restrict the X-axis from -5 to +5 ns. This eliminates any wrong pulses coming from electrical noise etc. Please note the "Mean" of the histogram which is 1.484 ns in the example below and the number of events "N" which is 168 in the example. You will use your values to do the following measurements.



# 3.2 Measurement 1: Speed of Cosmic Muons

To measure the speed of cosmic muons, the time of the signal of the lower plate can be subtracted with the time of the signal on the upper plate. This is effectively done in the histogram above and can be read under "Mean". Unfortunately, the cables do not have exactly the same length, which adds to the time measurement and causes a large uncertainty. To cancel this errors, two measurements must be made, one with a normal orientation of the detector, and one with the "upside-down" orientation, see below:



Normal orientation

Inverted orientation

For each measurement, take at least 200 events and note the "Mean" of each histogram. This takes about 15 minutes. To avoid any noise events far at left or right, restrict the X-axis to -5 ns to +5 ns by dragging the handles on the left and right of the histogram towards the middle.

After you wrote down the mean of the histogram, rotate the detector upside-down and clear the histogram bit hitting "CLR". Now repeat the measurement with 200 events and note the new mean. Then calculate the speed of the muons as a percentage of the speed of light ( $3 \times 10^8$  m/s), assuming that the detector plates are 10 cm apart:



## 3.3 Measurement 2: Direction of Cosmic Muons

Measure the direction of cosmic muons. The coincidence between the two scintillator plates effectively selects muons travelling through both plates, so it defines a direction. Let's rotate the scintillator plates in steps of 10 degrees and take events for 5 minutes each. Write down the number of events  $N_i$  for each measurement (i=0, 10, 20, ... 90).

## 3.4 Data Analysis

Analyze your data using Jupyter Notebook and Python. To create a new notbook, go to

#### https://jupyter.org

and click on "Try in your browser", then on JupyterLite, then on Notebook - Python:



You created your first notebook. Rename it (right mouse-click on the title) to something like For the first exercise "Speed-CosmicMuons.jpynb"

For the second exercise, create another notebook and name it similarly, for example, "Cosmics Angle.jpynb".

The notebook can contain two types of text, either plain text using the

Markdown language https://en.wikipedia.org/wiki/Markdown

or

Python code <a href="https://en.wikipedia.org/wiki/Python">https://en.wikipedia.org/wiki/Python</a> (programming language)

The cell type can be chosen by either selecting "Markdown" or "Code" on the top right selection box. Familiarize yourself with the notebook. Add some text and a bit of python code.

Click the right arrow ► to interpret the Markdown or execute the Python code.

#### 3.4.1.1 Measurement 1: Speed of Cosmic Muons

#### The Markdown code for the following text is:

```
# Measurement 1: Speed of Cosmic Muons
```

 $\star$  Take at least 200 events and note the "Mean" of each histogram (about 15 minutes).

\* Input mean time differences (in ns) for normal and inverted orientations

- \* Convert time from ns to seconds
- \* Calculate time difference
- \* Calculate muon speed
- \* Calculate speed as a percentage of light speed



As the next step, write some Python code that can calculate the speed of muons and the speed as a percentage of light speed.

You can ask ChatGPT to write some code for you. Go to <u>https://chatgpt.com</u> and make yourself a free account. Then try the following prompt:

Write an easy Python script that calculates the speed of muons based on user input of mean time differences for both normal and inverted orientations in nanoseconds.

The script should: 1. Convert these times to seconds. 2. Use a fixed distance of 0.1 meters between plates. 3. Calculate and display the delta time, speed of muons, and speed as a percentage of the speed of light. 4. Display the results to 2 decimal places.

Review the code generated by ChatGPT and adjust it as needed. Never trust the output of a Large Language Model, since there is no guarantee that it is correct. **Always double check!** 

#### 3.4.1.2 Measurement 2: Direction of Cosmic Muons

The Markdown code for the following text is:

```
# Cosmic Angle Measurement
```

\* Take data at angles 0, 10, 20, ... 90 degrees

Time-of-Flight Data Acquisition Lab

```
* Integrate counts for 5 minutes* Write down time / value pairs in notebook
```

You can execute Python code directly in the Jupyter Notebook, if you change the cell type from "Markdown" to "Code". For example "print(10 \* 5)" does the actual calculation 10 \* 5 and prints the result.



Now write some Python code to plot your array of ten X / Y values. If you are a Python expert, proceed immediately. If you don't know Python, ask ChatGPT to write some code for you.

Try the following prompt:

Write me some python code to plot an array of ten X/Y pairs. The X values represent an angle going from 0 to 90 degrees in steps of 10 degrees, the Y values represent count rates. Put ten discrete count rates from 0 to 100 into an array and normalize it from 0 to 1.

Review the Python Code from ChatGPT. Adjust titles, font sizes and colors if necessary and put your real count rates measured above into the code. Then run it.

The result should look like that:



Overlay the normalized counts with a  $cos(\theta)$  curve and  $cos^2(\theta)$  curve to see which one fits better. Figure out yourself how to get the Python code for that. As the result, you should get a plot similarly to this one:



If you cannot access ChatGPT, you can use this Python code directly and modify the <code>count\_rates</code> values:

```
import numpy as np
import matplotlib.pyplot as plt
# Define the X values representing angles from 0 to 90 degrees
x_values = np.arange(0, 91, 10)
x radians = np.deg2rad(x values) # Convert degrees to radians for trigonometric functions
# Define the Y values representing count rates from 0 to 100
count_rates = np.array([100, 98, 90, 79, 68, 56, 31, 24, 13, 17])
# Normalize count rates to range 0 to 1
normalized count rates = count rates / count rates.max()
# Calculate cosine and cosine squared values
cos_values = np.cos(x_radians)
cos squared values = cos values**2
# Plotting the data
plt.figure(figsize=(10, 6))
plt.plot(x_values, normalized_count_rates, marker='o', label='Normalized Count Rates')
plt.plot(x_values, cos_values, marker='x', linestyle='--', label='cos(\theta)')
plt.plot(x_values, cos_squared_values, marker='s', linestyle=':', label='cos^2(\theta)')
plt.title('Normalized Count Rates vs Angle with Cosine Overlays')
plt.xlabel('Angle (degrees)')
plt.ylabel('Normalized Values')
plt.legend()
plt.grid(True)
plt.show()
```

## 3.5 Discuss your measurement

Critically question your data, and think about the following questions:

- Why does the rate not drop to zero at 90 degrees?
- Why are the points not on a smooth line?
- If you measure again, will you get exactly the same points?
- How could the experiment be improved?