


Installation and Performance of Jungfrau Detectors at Diamond Light Source

John Matheson, Danny Axford, Anna Bergamaschi, Maria Carulla, Nicholas Devenish, Noemi Frisina, Matthew Hales, Viktoria Hinger, Ken Jones, Vadym Kedych, Andy Marshall, Aldo Mozzanica, Eva Gimenez-Navarro, James O'Hea, Dominic Oram, Robin Owen, Adam Prescott, Bernd Schmitt, Adam Taylor, Graeme Winter, Gary Yendell



Jungfrau




Jungfrau pixel architecture

Jungfrau is a modular detector family, based on tiling 0.5Mpixel modules. A module uses 4x2 ASICs, bump bonded to a single sensor. Pixel size is 75µm square. Sensors are 320µm Si. Very large dynamic range, from single photon to 10⁴ photons/frame, is enabled by integrating the input charge, with adaptive gain provided by switching feedback capacitors. Continuous Frame rates up to 2.2kHz are available and pixel output is encoded by a 16 bit number, giving data rates near 20Gb/s per module. Power consumption is near 40W per module. The integrating approach requires cooling the sensors to below -10°C.

Services for 9M

- Standard European 3-phase ac mains
- Chiller 1 (sensors): Lauda VC5000, -15°C flow, 1.3kW capacity, 3.2 bar pump
- Chiller 2 (power supplies): SMC HEC006A5B, 20°C
- Readout boards (ROBs) are fan cooled
- Dry air purge for cold compartment
- Copper network for slow control
- Fibre network for data, 850nm MM, 9 x 40Gb/s



Jungfrau-9M

Interlocks for 9M

Potential failure modes

- Cooling failure - overheating of front ends
- Dry air failure - condensation on front ends
- Condensation on ROBs, if cooled but unpowered

Sensors for Interlocks

Installed as standard by PSI:

- Sensirion SHT-35 sensor for temperature & RH


Added in DLS installation:

- K-type thermocouple
- Honeywell HIH-4000-003 sensor for RH
- Belimo 22PF-1UC flow meter on cold flow

Jungfrau 1M at DLS

Jungfrau-1M evaluated in Lab and on Diamond Beamline I24. Examples:

- Darks behaviour
- Single photons using X-ray tube (Cu-Kα, 6keV) (1)
- Direct beam from X-ray tube → gain switching (2)
- Diffraction from known samples (e.g. insulin) → linearity, dynamic range (3)




Although a switching discontinuity exists in raw pixel data, it is small under sparse illumination and doesn't affect integrated diffraction spots significantly.

For more details: J.Matheson et al., Journal of Synchrotron Radiation: Volume 33, Issue 2 Pages: 331-343 (20th February 2026)


Acceptance tests for 9M

Jungfrau-9M tested on Diamond Beamline I24

- Power supply currents
- Compare noise maps before/after shipping (1)
- Check reconstruction of data streams using LaB6 rings (2)
- Single photons calibration check on Beamline I24 (12.4 keV) with glycerol as an amorphous scatterer (3)




Hardware Layout for 9M




- "Grace Hopper" - QuantaGrid S74G-2U server
 - [Hyperscale Products | QCT- Quanta Cloud Technology](#)
 - <https://www.nvidia.com/en-gb/data-center/grace-hopper-superchip/>
- Trigger from beamline motion using Quantum Detectors PandaBox

DAQ Software for 9M



- Single Grace Hopper Uses 57% of GPU capacity at 1kHz frame rate
 - Measured using the NVIDIA-SMI utility
- Software optimisation expected to allow 2kHz running
- Bottleneck is memory access for darks subtraction / gain correction

9M in use



Jungfrau-9M detector on Beamline I24 (1)

Section of a diffraction image from human deacetylase (1kHz) (2)

Resulting electron density obtained from a single crystal at 100 K (3)