

Doppler shift lifetime measurements using the TIGRESS Integrated Plunger at ISAC-II/TRIUMF



K. Starosta for TIGRESS/TIP collaboration

Department of Chemistry
Simon Fraser University

May 29, 2017




Selected goals of nuclear science research



- Understand the mechanisms of shell evolution in medium-mass and heavy nuclei as a function of isospin.
- Develop a theoretical framework that is able to make accurate predictions of nuclear properties.


Studying nuclear structure using the electromagnetic force

- The electromagnetic force provides a convenient non-intrusive probe of nuclear systems bound by the strong force.
- Lifetime measurements using gamma-ray spectroscopy provide:
 1. An observable sensitive to nuclear structure.
 2. A useful benchmark for nuclear model calculations.




$$\begin{aligned}\tau(E2; 2_1^+ \rightarrow 0_1^+) &= \lambda(E2; 2_1^+ \rightarrow 0_1^+)^{-1} \\ \lambda(E2; 2_1^+ \rightarrow 0_1^+) &\propto E(2_1^+)^5 \times B(E2; 2_1^+ \rightarrow 0_1^+) \\ B(E2; 2_1^+ \rightarrow 0_1^+) &= \frac{1}{5} \langle 2_1^+ || E2 || 0_1^+ \rangle^2 \propto \beta^2\end{aligned}$$

The Recoil Distance Method with TIP and TIGRESS




The TIGRESS Integrated Plunger (TIP) device



P. Voss et al. Nucl. Inst. and Meth. A 746 (2014) 87, P. Voss et al. Phys. Proc. 66 (2015) 524.

TIP with the CsI(Tl) wall




TIP commissioning experiment details

- Measurement of the $2_1^+ \rightarrow 0_1^+$ lifetime in ^{84}Kr .
- Previous Coulex experiment reports $\tau = 5.84 \pm 0.18$ ps.
- TIP and 11 TIGRESS detectors in a 3/5/3 configuration with 24-element CsI(Tl) wall for particle identification (PID).
- Excited state populated via partially unsafe Coulex reaction (beam energy of 250 MeV or 2.979 MeV/u).
- A total of 13 target/degrader separation distances from 20–400 μm were analyzed.
- Data analysis via a comparison to Geant4-simulated lineshapes developed for low-statistics experiment analysis.

T. J. Mertzimekis et al. Phys. Rev. C 64 (2001) 024314.


Geant4 simulation framework



- Coulomb excitation followed by gamma-ray decay.
- Analytic solutions for single step $E2$ process (Coulex kinematics, angular distributions, etc.) with track weighting to handle thick target integration.
- Gamma-ray sensitive detectors ported from GRIFFIN/TIGRESS code originating from Guelph.

Adrich et al. Nucl. Inst. and Meth. A 598 (2009) 454, Alder et al. Rev. Mod. Phys. 28 (1956) 432.

Geant4-facilitated data analysis: Doppler-shift factors



Histogram analysis using the likelihood ratio χ^2

- For a Poisson likelihood function, the likelihood ratio χ^2 is given by

$$\chi^2 = 2 \sum_{i=1}^k y_i - n_i + n_i \ln(n_i/y_i)$$

where y_i is the model and n_i the observed data.

- Pros:
 - Versatile (goodness of fit, point estimation, error analysis).
 - Control over parent distribution.
 - Minimizing likelihood ratio $\chi^2 \equiv$ maximizing likelihood function.
 - **No variance estimation!**
- Cons:
 - Non-linear.
 - ...

Data analysis procedure


- For a given input lifetime τ , simulate gamma-ray spectra grouped by the Doppler-shift factor at all distances.
- Model data using

$$y_i = \alpha_0 s_i + \alpha_1 + \alpha_2 \operatorname{erfc} \left(\frac{i - c}{w\sqrt{2}} \right),$$

where s_i is the simulated data and the α 's are free parameters.

- Minimize $\chi_{d,g}^2$ for each distance and group.
- Minimum in total $\chi^2 = \sum_{\text{dist.}} \sum_{\text{gr.}} \chi_{d,g}^2$ corresponds to best fit lifetime τ_{\min} .

Best fit lifetime from χ^2 analysis



Comparison of fitting methods

1. Likelihood ratio χ^2


$$\chi_{\lambda}^2 = 2 \sum_{i=1}^k y_i - n_i + n_i \ln(n_i/y_i)$$

2. Neyman χ^2


$$\chi_N^2 = \sum_{i=1}^k \frac{(n_i - y_i)^2}{n_i}.$$

- Probably the most common fitting statistic (ROOT default for histogram fits, for example).
- Uses data to estimate variance.
- Assumes bin error can be estimated by $\sqrt{n_i}$: what if $n_i = 0$?


Lineshape fits for ^{84}Kr , 60 μm : all data



Lineshape fits for ^{84}Kr , 60 μm : 1% of data



Simulated lineshapes: groups 1 and 5




Best fit lifetime: 5.880 ± 0.013 (stat.) ± 0.070 (sys.) ps

Literature value: 5.84 ± 0.18 ps [Mertzimekis 2001]

Commissioning experiment summary

- Systematic uncertainties from the following 3 sources were identified:
 1. Transitions from higher-lying (feeding) states,
 2. Misalignment of the target and degrader foils,
 3. Choice of fit range for the χ^2 analysis.
- No deorientation effect observed in the data.
- Final reported lifetime: $\tau = 5.880 \pm 0.008$ (stat.) ± 0.070 (sys.) ps.
- Excellent agreement with literature value of 5.84 ± 0.18 ps with factor of ~ 2 reduction in uncertainty.
- A robust and flexible framework has been developed for the planning and analysis of RDM experiments using TIP.
- Paper for submission to Nucl. Inst. and Meth. A in preparation.

TIP RIB experiment motivation




TIP RIB experiment details

- Performed December 9–14, 2015.
- Measurement of the $2_1^+ \rightarrow 0_1^+$ lifetime in ^{94}Sr .
- Previous fast timing experiment reports 10 ± 4 ps.
- TIP and 16 TIGRESS detectors with 24-element CsI(Tl) wall.
- Running conditions same as commissioning experiment.
- Except the beam rate: factor of $\sim 10^4$ lower for RIB.
- Data was recorded at three target/degrader separation distances: 50, 100, 150 μm .


Mach et al. Nuc. Phys. A 523 (1991) 197.

Particle identification with RIBs




(a) Electrons from β^- decay (b) Al recoils from Coulex reaction


Best fit lifetime from χ^2 analysis




Simulated lineshapes: groups 1 and 5



$B(E2 : 2_1^+ \rightarrow 0_1^+)$ measurement



Impact of high-precision $B(E2 : 2_1^+ \rightarrow 0_1^+)$ measurement



Summary of results

Current work

Lifetime τ $7.80_{-0.4}^{+0.5}$ (stat.) ± 0.07 (sys.) ps

$B(E2; 2_1^+ \rightarrow 0_1^+)$ $0.0254_{-0.0014}^{+0.0015}$ (stat.) ± 0.0002 (sys.) e^2b^2

Mach (1991)

Lifetime τ 10 ± 4 ps

$B(E2; 2_1^+ \rightarrow 0_1^+)$ $0.020 \pm 0.008 e^2b^2$

- A robust and flexible framework has been developed for the planning and analysis of RDM experiments using TIP.
- Details of simulation framework, data analysis, and results for the ^{84}Kr commissioning experiment for submission to Nucl. Inst. and Meth. A in preparation.
- ^{94}Sr lifetime measurement paper submitted to Phys. Rev. C.

Acknowledgments

TIP Design

R. Henderson, TRIUMF

Simon Fraser University

C. Andreoiu, D. Cross, T. Domingo, F. H. Garcia, J. L. Pore, U. Rizwan, K. Starosta, P. Voss, J. Williams

SFU Science Machine & Electronics Shops

R. Holland, P. Kowalski, J. Shoults, K. Van Wieren

TRIUMF

G. Ball, N. Bernier, R. Caballero-Folch, S. Cruz, L. Evitts, A. B. Garnsworthy, G. Hackman, S. Hallam, J. Henderson, R. Krücken, E. MacConnachie, M. Moukaddam, O. Paetkau, C. Pearson, P. Ruotsalainen, J. Smallcombe, J. K. Smith, M. Williams

University of Guelph

B. Olaizola, E. T. Rand, C. E. Svensson

Universidad Nacional Autónoma de México

E. Padilla-Rodal,

University of Tokyo

K. Wimmer

University of Toronto

T. E. Drake

St. Mary's University

R. A. E. Austin



Theory references

- Kumar 1985: Quadrupole plus pairing Hamiltonian.
- Skalski 1997: Nilsson-Strutinsky method with the Woods-Saxon potential.
- Lalazissis 1995: Relativistic mean field.
- Buchinger 1994 (a): Finite range liquid drop.
- Buchinger 1994 (b): Strutinsky energy theorem.

A. Kumar and M. R. Gunye, Phys. Rev. C 32, 2116 (1985).

J. Skalski, S. Mizutori, and W. Nazarewicz, Nuclear Physics A 617, 282 (1997).

G. Lalazissis and M. Sharma, Nuclear Physics A 586, 201 (1995).

F. Buchinger *et al.*, Phys. Rev. C 49, 1402 (1994).

TIP commissioning experiment details

^{84}Kr properties

E_γ	881.615 keV
$\tau_{\text{lit.}}$	5.84 ± 0.18 ps

Plunger setup.

	Material	Thickness [mg/cm^2]	Thickness [μm]
Target	Al	1.07 ± 0.04	3.96 ± 0.16
Degrader	Cu	3.90 ± 0.16	4.35 ± 0.18

Beam properties

Beam energy	250 MeV
Safe Coulex	200 MeV
Rate	$\sim 2 \times 10^8$ pps

TIP RIB experiment details

^{94}Sr properties

E_γ 836.9 keV

$\tau_{\text{lit.}}$ 10 ± 4 ps

Plunger setup.

	Material	Thickness [mg/cm^2]	Thickness [μm]
Target	Al	1.09 ± 0.04	4.05 ± 0.17
Degrader	Cu	3.69 ± 0.15	4.12 ± 0.17


Beam properties

Beam energy 280 MeV


Safe Coulex 227 MeV

Rate $\sim 2\text{--}5 \times 10^4$ pps


Thick target integration with Geant4




Feeding analysis



$$\begin{aligned}
 N_5(t) &= N_4(0) \left(1 - e^{-\lambda_4 t}\right) \\
 &+ \left[\sum_{i=1}^3 N_i(0) \frac{\lambda_i}{\lambda_4 - \lambda_i} \right] \left(e^{-\lambda_4 t} - 1\right) \\
 &+ \sum_{i=1}^3 \left[N_i(0) \frac{\lambda_4}{\lambda_4 - \lambda_i} \left(1 - e^{-\lambda_i t}\right) \right]
 \end{aligned}$$



Deorientation effect



Analysis of low-statistics data sets

Challenges associated with the comparison of simulated data sets from Geant4 to low statistics data:

- Typical least squares analysis requires an estimate of the variance either from the model (Pearson χ^2) or from the data (Neyman χ^2).
- These are derived under the implicit assumption of a Gaussian error distribution.
- With few counts, errors are poorly estimated.
- Even worse, if there are 0 counts in the data, the very commonly used Neyman χ^2 statistic is undefined!
- Inconsistent normalization (which depends on choice of statistic!) without inclusion of explicit normalization parameter.

Basic quantities for histogram analysis

Consider a histogram with k bins labeled by the index i running from 1 to k and a model with J parameters labeled by index j . Define the following quantities:

- $n_i =$ number of events in bin i .
- $\mathbf{n} = (n_1, n_2, \dots, n_k)$.
- $\boldsymbol{\alpha} = (\alpha_1, \alpha_2, \dots, \alpha_j)$.
- $y_i =$ number of events predicted by the model (via $\boldsymbol{\alpha}$) in bin i .
- $\mathbf{y} = (y_1, y_2, \dots, y_k)$.

Baker and Cousins. Nucl. Inst. and Meth. A 221 (1984) 437.

Likelihood ratio χ^2 derivation

Assume the data in each bin are independently Poisson distributed. Then the likelihood function \mathcal{L} is given by

$$\mathcal{L}(\mathbf{y}; \mathbf{n}) = \prod_{i=1}^k \frac{y_i^{n_i} \exp(-y_i)}{n_i!}.$$

Let \mathbf{m} be the true (unknown) values of \mathbf{n} . The likelihood ratio Λ is

$$\Lambda = \frac{\mathcal{L}(\mathbf{y}; \mathbf{n})}{\mathcal{L}(\mathbf{m}; \mathbf{n})},$$

and the likelihood ratio test theorem says that the “likelihood chi-square”

$$\chi_{\Lambda}^2 = -2 \ln \Lambda = -2 \ln \mathcal{L}(\mathbf{y}; \mathbf{n}) + 2 \ln \mathcal{L}(\mathbf{m}; \mathbf{n})$$

asymptotically obeys a χ^2 distribution.

Baker and Cousins. Nucl. Inst. and Meth. A 221 (1984) 437, Wilks. Ann. Math. Stat. 9 (1938).

Likelihood ratio χ^2 derivation

Replace the unknown \mathbf{m} with its bin-by-bin maximum likelihood estimator \mathbf{n} . A bit of algebra yields

$$\begin{aligned}\chi_{\lambda}^2 &= -2 \ln \mathcal{L}(\mathbf{y}; \mathbf{n}) + 2 \ln \mathcal{L}(\mathbf{m}; \mathbf{n}) \\ &= -2 \ln \mathcal{L}(\mathbf{y}; \mathbf{n}) + 2 \ln \mathcal{L}(\mathbf{n}; \mathbf{n}) \\ &= 2 \sum_{i=1}^k y_i - n_i + n_i \ln(n_i/y_i).\end{aligned}$$

Some attractive features:

- It's a χ^2 statistic (familiarity, versatility, error analysis).
- No variance estimation.
- Self-normalizing.
- A clear way to handle bins with $n_i = 0$.
- Minimizing χ_{λ}^2 is equivalent to maximizing $\mathcal{L}(\mathbf{y}; \mathbf{n})$.