

Positron emission particle tracking using a modular positron camera

*D J Parker, T W Leadbeater, X Fan,
M N Hausard, A Ingram and Z Yang*



UNIVERSITY OF
BIRMINGHAM


Positron Imaging Centre

Positron emission tomography (PET):

Mapping concentration of radioactively-labelled fluid

PET scanner consists of rings of many small detectors, operating in coincidence to detect the pairs of back-to-back γ -rays from positron annihilation.

After detecting millions of such events a 3D tomographic map of tracer concentration can be reconstructed




Engineering PET: Gamma rays are penetrating – can observe labelled fluid inside process vessels

But PET is slow – requires detection of ~ million coincidence pairs


For fast dynamic information use **Positron emission particle tracking (PEPT):**

Introduce a single labelled particle, and locate it frequently




Detection

Detection of gamma rays using two large position sensitive detectors.



Reconstruction

Two rays detected in coincidence define line along which particle lies.




Particle Location

After several events, tracer can be located via triangulation.

Currently labelling tracer particles down to 100µm diameter

Can locate tracer particle to within 1mm every 1ms

PSD1: Original Birmingham Positron Camera




- Originally designed to image lubricant in operating jet engine
- Operational 1984
- Pair of gas filled MWPCs, sensitive area $600 \times 300 \text{mm}^2$

- Cathode planes have $50 \mu\text{m}$ lead strips
- Delay line readout
- Each chamber contains a stack of 20 such assemblies
- Total quantum efficiency for 511keV photons 7%

- Spatial resolution 8mm FWHM (+long tails)

Useful count rate limited to around 3000 events/s due to
dead time in readout
random coincidences (resolving time 12.5ns)



PSD5: “New” Birmingham Positron Camera




- Installed summer 1999
- Commercially available medical system (ADAC Forte)
- Cost £0.3M
- Pair of gamma camera heads on rotating gantry; separation 250-800mm

Each head contains NaI(Tl) crystal 590x470x16mm³,
55 PMTs each connected to its own ADC;
single board computer




	Old camera	New camera
Spatial resolution	8mm	4-6mm
Singles sensitivity	7%	23% total/16% photopeak
Energy resolution	None	15% FWHM
Resolving time	12.5ns	7.5ns
Max true coincidence rate	5kcps	100kcps

PEPT data from spouted fluidised bed



A typical particle track
This track represents
five seconds worth of
particle motion within a
spouted fluidised bed.

The amount of
time spent by the
particle within a
particular region of
the bed, or
occupancy, can be
related to the
density of the bulk.



Analysis of the
particle track allows
production of vector
plots.

Examples of standard PEPT



Scraped-surface heat exchanger



Fluidised bed

PSD7: Third generation Birmingham systems developed from redundant medical PET scanners

Since 2002 PIC has acquired
4 complete PET scanners
and components from two others


These are inherently modular, and can in principle
be reconfigured in different geometries.

Some are being used/developed for PET
(uniform rings)
and others for PEPT



They are based on segmented BGO blocks ($60 \times 50 \times 30 \text{mm}^3$) cut into 8×4 or 8×8 arrays, and viewed by 4 PMTs

These are grouped in modules (“buckets”) comprising 4 blocks plus associated electronics




For PEPT, use these modules (in appropriate geometry) and existing coincidence processors, etc

Just need to store data in list mode


Transportable Modular Positron Camera for PEPT

16 detector modules (extendable to 32)

Can be configured in different geometries:



Count rates (close geometry, central source)




Curves are (absolute) predictions of simple model, assuming $0.75\mu\text{s}$ deadtime in each bucket, and $1\mu\text{s}$ deadtime in data storage system

- Latest data acquisition system (T Leadbeater poster) does much better

Effect of gaps on sensitivity

Example: two rectangular arrays, 2 boxes wide x 4 boxes high, separated by 500mm
overall dimensions of each array 666 x 345 mm² (comparable to ADAC)
but 60% of area is “gap”



Gaps have drastic effect on coincidence sensitivity


But quantum efficiency of 30mm BGO is ~3x higher than 16mm NaI(Tl)

So overall, average sensitivity of modular camera is still higher


On-plant PEPT study : BP, Hull

(240km from Birmingham)


750mm diameter fluidised bed, with central dividing baffle + different air supplies each side of baffle




4 banks of detectors (detector separation 1.2m) give FOV shown





Analysis in terms of movement between four quadrants



With asymmetric gas flow to promote anticlock circulation



Casting of liquid aluminium: PEPT tracking of small alumina inclusion



Summary

Modular camera enables PEPT studies on large systems and outside the lab

Gaps between blocks seriously reduce sensitivity

Count rate was limited by data storage

– improved system (Leadbeater) gives $>1 \text{ M s}^{-1}$

Acknowledgements

This work was made possible thanks to financial support from EPSRC and Yorkshire Forward