



UNIVERSITÀ
DEGLI STUDI
DI TORINO



LGAD-based detectors for monitoring therapeutic proton beams

A Vignati, M Abujami, G Borghi, F Ficarella*, C Galeone, S Garbolino,
S Giordanengo, O Hammad Ali*, OA Martì Villarreal, F Mas Milian, G Mazza,
RJ Wheadon, R Cirio, R Sacchi, V Monaco*

12th International Conference on
**POSITION SENSITIVE
DETECTORS**



**Fondazione Bruno Kessler (FBK)*



Hosted by
**UNIVERSITY OF
BIRMINGHAM**



Acknowledgments



This work was supported by the INFN Gruppo V (MoVe-IT project), by the European Union's Horizon 2020 Research and Innovation funding program (Grant Agreement no. 669529 - ERC UFSD669529), by the Ministero della Ricerca, Italia, PRIN 2017 (project 2017L2XKTJ – 4DinSiDe). We acknowledge the proactive collaboration of the Fondazione Bruno Kessler to this project. We also thank the CNAO Foundation (Pavia) and the Centro di Protonterapia (Trento) for the beam time accorded.




Particle Therapy



dose released to the tissues




Dose Delivery System




CNAO – Pavia IT

Dose Delivery System



A correct dose distribution needs a precise **monitoring** of the beam parameters:
fluence, position, shape, energy




Modelling and Verification for Ion beam Treatment planning

Development of two prototypes of innovative detectors based on silicon (UFSDs) for beam monitoring:

1. to **directly count** individual protons
2. to measure the **beam energy** with time-of-flight techniques

Ultra Fast Silicon Detectors (UFSD)



Based on the Low Gain Avalanche Diode (LGAD) technology

- Thin p⁺ gain layer implanted under the n⁺⁺ cathode
 - Controlled low gain ($\sim 10\text{-}30$)
 - Gain increases with bias voltage
 - High Signal/Noise ratio
- Small thickness of active volume 50 μm
 - Reduced beam perturbation, signal steepness
 - Loss of signal due to reduced thickness compensated by the internal gain



*N Cartiglia
(Talk on
Tuesday)*

Short signal duration ($< 2 \text{ ns}$)

→ **particle counting**

Excellent time resolution (tens of ps)

→ beam **energy measurement** from time-of-flight





Synchrotron

FWHM 10 mm

Flux $10^8 - 10^9$ p/s in spills

Intensity 20 - 50 - 100%

Energy 62 - 227 MeV



Cyclotron

FWHM 3 - 7 mm





Flux $10^6 - 10^{10}$ p/s

Current 1 - 320 nA


Energy 68 - 228 MeV

Particle Counting










Dedicated fast ASIC chip (ABACUS)



- Dead time → ~ 10 ns
- Counting accuracy at 1%
- 110 nm CMOS technology
- Chip area = 2 × 5 mm²
- 24 channels
- CSA dynamic range: 4 - 150 fC



Dedicated fast ASIC chip (ABACUS)





Dead time $\rightarrow \sim 10$ ns

FPGA


- FE initialization
- Pulse counting
- Pileup correction

CSA dynamic range: 4 - 150 fC


Threshold scans with 4 different proton beam energies (70 MeV, 125 MeV, 179 MeV, 228 MeV)



Signal amplitude distributions - channel 7



MPV vs energy - channel 7



$$E_{beam} \propto \beta^2$$


Tests with the ABACUS testboard at




TIFPA



Signal amplitude distributions - channel 7




MPV vs energy - channel 7



60 MeV – 230 MeV
30 mm – 320 mm



$$E_{beam} \propto \beta^2$$



Particle Counting


Pile-up inefficiency

Mitigation of saturation effects




Particle rate estimated using the charge measured with a PTW pin-hole ionization chamber and assuming a paralyzable saturation model.

Correction method based on the correlation of logical pulses from two neighbouring strips (paralyzable model of saturation effects assumed)




Final Counting Sensor Prototypes

- Strip segmentation (strip area $\sim 3 \text{ mm}^2$)
- Area $2.7 \times 2.7 \text{ cm}^2$
- 146 strips (144 with gain, 2 no gain)



OA Marti Villarreal
(Talk on Tuesday)





Energy Measurement





Energy Measurement Telescope

- Determine the energy through measurements of **time-of-flight** of coincident protons
- Required tolerance on the range uncertainty in water **< 1mm**, i.e. precision **< 4 ps** needed at 230 MeV at 1 m distance





Timing Prototypes and Readout





11 strips,
pitch 590 μm
50 μm active thickness
(thinned down to a total
thickness of 70 μm)



Dedicated
readout board
for time
measurement



Time-of-Flight (TOF) measurement - Coincidences




Time-of-Flight (TOF) measurement - Coincidences



constant fraction algorithm to reduce the time walk effect

t_1

threshold



$\Delta t = t_1 - t_2$


t_2

Distribution of Δt

- Double gaussian fit to identify the peak of true coincidences from the combinatorial background


CNAO

181.68 MeV @ 400 mm



TIFPA PROTONTERAPIA
TRENTO

182.8 MeV @ 400 mm




Distribution of Δt

- Double gaussian fit to identify the peak of true coincidences from the combinatorial background


CNAO

181.68 MeV @ 400 mm





TIFPA PROTONTERAPIA
TRENTO

182.8 MeV @ 400 mm



Energy measurement from TOF


→ Interpolation of the measured time difference as a function of the distance with K_0 and *time offset* as parameters




A **calibration** is needed to remove the systematic errors due to the experimental setup

Preliminary Results

CNAO



PROTONTERAPIA TRENTO



UFSD are a promising technology for beam qualification and monitoring in Particle Therapy

- Excellent time resolution: **real-time measurement of the beam energy**
- Short signal duration: **single particle counting**

➔ *Patent for energy measurement prototype*

➔ *Final readout electronics for counting prototype*

➔ *Radiation hardness*

➔ Radiobiological experiments

➔ **Aiming at clinics!**



Thank you

